

Operation Light






KD-2036L-3



Mail: contacto@healthcorp.ec
 Celular: 0999230608

Mobile Operation Light



Features:

- ♦ AC supply,DC to choose
- ♦ Long arm,wide scope of work
- ♦ LED source lighting,no temperature rise at work scope

KD-2036L-3

Rated voltage	Bulb power	Illumination	Color T	Lamp life	Weight	Package
AC90~240V 50Hz~60Hz	LED3Wx36	> 65000Lux	4500±500K	50,000H	33.5Kg	1140×625×316mm