

LIPOZERO

SPECIAL BODY & FACE LINE CARE

LIPOZERO

SPECIAL BODY & FACE LINE CARE

- Improvement of body-shape ✓
- Improvement of facial-line ✓
- Improvement of obesity ✓
- Improvement of cellulite ✓
- Regeneration of collagen ✓

SPECIFICATION

Technologies	RF, Ultrasound Cavitation, Low Level Laser, LED, Vacuum
Power	AC100~240V, 50/60Hz
Weight	8.5kg
Dimension(WxDxH)	371 x 414 x 241(mm)



Healthcorp
MEDICAL EQUIPMENT

Mail: contacto@healthcorp.ec
Celular: 0999230608

Eunsung




LIPOZERO

SPECIAL BODY & FACE LINE CARE

LIPOZERO is the solution for total body-line management!
Two types of handpieces can treat from the face and body at the same time!

Body & Face Line Smart Care




The combination of a different technologies for improved fat reduction!

Reduces the fat and cellulite at the same time, leading to improvement of over all healthier body

Solves out the body care & Anti-aging at the same time

Left Handpiece



1 Synergy effect through 4 different technologies

Able to perform RF, Cavitation, Low level laser, and LED either all together or separately.

2 Deduction of fat cells

By liquefying the fat and let them easily burned thorough physical activity.

3 Acceleration of metabolism

Able to excrete the fat cells in a short time by revitalizing metabolism of the inner body.

4 Systematic care of individual body parts

Able to concentrate treatment on specific body parts.

Right Handpiece

1 Synergy effects through 3 different technologies

RF, Vacuum, LED can be performed either together or separately.

2 Able to perform on both body and face


Having diverse size of handpiece that enable treatment on both body and face.

3 Acceleration of lipolysis and metabolism

Deep heat from RF energy induce lipolysis and accelerate metabolism.

4 Stimulation of cellulite and massage effect

Able to physically stimulate cellulite using strong vacuum while stimulating muscle tissue for massaging effect at the same time.



Five functions of LIPOZERO



RF

Multi polar RF energy with 0.8Mhz generate deep heat that will draw many effects, such as deduction of subcutaneous fat, acceleration of metabolism, collagen regeneration and so on.



Vacuum

By utilization of vacuum, RF output can effect down to the subcutaneous tissues.



Cavitation

Micro-cavities (bubbles) generated by low frequency ultrasound will implode and destroy fat cells.



Low Level Laser

Low Level Laser with 660nm, creates temporary pore on the fat cells and leads to the excretion of triglycerides from the tissues.



LED

Using Red & Blue LED separately or together, able to treat diverse skin problem.